

12/1996

BioGlo™ Sterile Fluorescein Strips

As per the manufacturer, this finished product is subject to the labeling and material information provisions of the Federal Food, Drug and Cosmetic Act. Therefore, it is their belief that a Material Safety Data Sheet for this product is not required under the OSHA Hazard Communication Standard.

888-EYE DROP (393-3767) • Ph: 001-909-476-8394 • Fax: 001-909-948-2193 • Email: orders@hubrx.com



9339 Charles Smith Ave., Bldg. 150 • Rancho Cucamonga, CA 91730

July 2, 1999

The ingredients used in the BioGlo™ Fluorescein Strips are:

1. Distilled water q.s
2. Fluorescein Sodium U.S.P.
3. Sodium Carbonate IP
4. Whatman Filter Paper grade #1

If I can be of any further service, please don't hesitate to contact me.

Jerry Rosenstein, O.D.
President
HUB Pharmaceuticals